




THE HADA

Crosslinked
Hyaluronic Acid



THE HADA Crosslinked Hyaluronic Acid

HA core Technology

Preserved Natural Entanglement Technology

The condensed story of the science behind PNET™

The unique and patented stabilization technique ensures that the natural cross-linked and entangled HA network is kept in place by introducing a limited number of synthetic cross-links, i.e. minimal modification.



THE HADA Crosslinked Hyaluronic Acid



THE HADA

Crosslinked
Hyaluronic Acid

-Medical grade HA (Crosslinked Hyaluronic Acid)

-HA volume : 50mL / vial

-HA concentration : 24mg / mL

-BDDE < 0.2ppm

-Gel type: Mono-phasic

-Package : 1vial / box



/ Safety

- Medical grade HA
- Minimal cross linker usage
- No side effects
- Minimal modification

/ Convenience

- Absorbable material : HA
- Easy to remove
- Soft gel type
- Short downtime

/ Economics

- Reasonable price
- Various use of HA
- Long lasting effect
- Outstanding lifting capacity

Advantages

Easy to use and produces natural results

Safe due to non-animal origin and low BDDE residue

Long-lasting effect by slow and uniform decomposition

Optimum results with low usage volume

