

# ELOHA *Plus*

Crosslinked Hyaluronic Acid  
with Lidocaine

HELOHA *Plus*



Crosslinked Hyaluronic Acid with Lidocaine

*Hello  
Hyaluronic Acid*





*ELOHA Plus is a new concept filler that shows high elastic modulus and excellent cohesion.*

# ELOHA *Plus*

Crosslinked Hyaluronic Acid  
with Lidocaine







# ELOHA<sup>Plus</sup>

Crosslinked Hyaluronic Acid  
with Lidocaine


Improving Wrinkle  
Volumizing  
Contouring  
Corrective Treatment for Asymmetry  
Improving Skin Quality





# ELOHA*Plus*

Crosslinked Hyaluronic Acid  
with Lidocaine






High Elasticity

High Cohesion

High Convenience

- Minimal Modification
- Outstanding Lifting Capacity
- High Viscosity
- Long Lasting Effect
- EC Certified

Name	Figures	Product Specification				Needle Specification		Property		
		HA	PHASIC (GEL TYPE)	BDDE	CONTENTS	GAUGE (UTW)	QUANTITY	VOLUME	LIFTING	COHESIVE
ELOHA Plus VOL.1		24mg/mL	MONO	<0.02ppm	1mL*1syringe	27G	2EA	<div><div></div></div>		
ELOHA Plus VOL.2		24mg/mL	MONO	<0.02ppm	1mL*1syringe	27G	2EA	<div><div></div><div></div></div>		
ELOHA Plus VOL.3		24mg/mL	MONO	<0.02ppm	1mL*1syringe	27G	2EA	<div><div></div><div></div><div></div></div>		

*ELOHA Plus is a new concept filler that shows high elastic modulus and excellent cohesion.*

# ELOHA *Plus*

Crosslinked Hyaluronic Acid  
with Lidocaine





# ELOHA<sup>Plus</sup>

Crosslinked Hyaluronic Acid  
with Lidocaine

## Cross Linking Process

Ideal cross linking: 100% of BDDE participate in cross linking, no not-linked BDDE exists.  
Some BDDE partially participated in cross linked(open cross linked), reported to cause adverse event.


USA FDA guide line: below 2ppm of open cross linked/non-cross linked BDDE.

## PNET(Preserved Natural Entanglement Technology)

Increase cross linking efficacy(Low BDDE) lead to minimize Open-cross linked BDDE.

Inactive open-cross linked BDDE through 1st cleaning process.

Eliminate free BDDE again in 2nd cleaning process.







# The Injection Area of

# ELOHA<sup>Plus</sup>

Crosslinked Hyaluronic Acid  
with Lidocaine

## ELOHA Plus Vol.1

- Forehead Lines
- Crow's Feet
- Neck Wrinkles
- Lower Eyeid Contouring



## ELOHA Plus Vol.2

- Nasolabial Folds
- Laugh Lines
- Forehead, Nose, Cheek
- Chin Augmentation



## ELOHA Plus Vol.3

- Nasolabial Folds
- Forehead, Nose, Cheek
- Chin Augmentation
- Glabellar Lines

