


iLLUMA PN Booster

Peptide ingredients and HA make
the skin elastic and healthy.
iLLUMA makes your skin bright and radiant.

iLLUMA PN Booster

Main Ingredients

- + PN Polynucleotides
- + Peptide
- + Hyaluronic Acid




iLLUMA
PN Booster

Add Skin Science to beauty

'iLLUMA PN Booster'
improves aging skin by strengthening
the skin barrier.

It is a skin booster that provides
moisture to the epidermis and dermis
and has an excellent effect on
repairing thin and damaged skin.




The unique and patented stabilization
technique ensures that the natural cross-
linked and entangled HA network is kept in
place by introducing a limited number of
synthetic cross-links, i.e. minimal
modification.

iLLUMA PN Booster is
an innovative method to fight
all age-related skin changes.

iLLUMA PN Booster

The PN polynucleotides (DNA + RNA) included
in the composition have a wider spectrum of
action compared to regular PDRN:

- Activate the general metabolism of the cell
- Normalize local skin immunity
- Stimulate the production of structural components of the dermis
- Block photo-damage of their own cell DNA, reducing the number of atypical cells
- Regulate melanogenesis by stabilizing cell membranes, which significantly reduces the sensitivity of cells to the damaging effects of chemicals and radiation.



The iLLUMA PN Booster also contains 2 anti-aging peptides:

- Acetyl hexapeptide-8 is specifically designed to combat facial wrinkles. It effectively reduces wrinkles on the neck and around the eyes. The peptide has proven to not only effectively reduce existing wrinkles, but to also prevent the appearance of new ones.
- Palmitoyl pentapeptide - 4 signal peptide that stimulates the restoration of the structural elements of the dermis - collagen, elastin, fibronectin, glycosaminoglycans.

/ Features

- Skin repair & hydration
- Facial wrinkle removal
- Strengthening of barrier properties and strengthening of local immunity of the cells
- Reduction of pigmentation of various etiologies and normalization of melanogenesis
- Stimulation of collagen and elastin synthesis

/ Interval

- 1 treatment per month,
between 1 to 3 sessions in total.



iLLUMA PN Booster